

Brochure for the Lt - light-thermo-therapy hand-appliance
Light-thermo-therapy method after Dr. Med. APUL

Content:

- At - operation
- B - application-areas
- C - contraindications
- D - piece of information and precautions
- E - side-effects
- F - interactions with other drugs
- G - type of the application and dosage
- H - Überdosierung and measures
- I - storage
- J - presentation-form and package-content
 - Manufacturers and distribution
 - Authorization-date and patent-number EP 1005 880
 - DP 19856002 - CE 0494 Nichtinvasive committees medical class gets 1.

.At - composition operation:

- TL 12° NGES as LAFER, with ca.700 facet-mirrors,
Light-sources 1 cm removals 170 mWs per 1 cm²s. Deliver
- Light-radiation leaves the appliance, that with improper application to combustion and Fire-generation can lead.
- , In order to hand over the light-energy at physical cells, , applies, this appliance was specifically produced.
- Our life-section is connected with many Stress; about the energy-loss interconnected with Stress to supplements, is this one treatment, that very much effective is.
- The treatment works circulation-promoting.
- The unspzifsche Muskuloskeleteal backaches are removed application.
 - Decrease the Muskelverspannungen. The increase of the pains can occur at the start.
- Pubs or regional Hyperthermie application.
- The application is effective also with depressions.

B - application-areas:

- To the prevention with tiredness, Stress, performance-waste, ages troubles,
 - and increase the vitality and the body-own immune-system.
 - saisonalabhängige depression (=SAD), easy until heavy-exactly-y depressive episode of all type,
- Fear-conditions, panic-stricken disturbances. Fatigue syndrome.
Chronic tiredness syndrome = CMS.
- Chronic rheumatic troubles forms circle.
 - Hals-WS, Brust-WSs, loins - spinal columns syndrome.
 - and Spondylitis ankylosans, Arthrose all type,
 - the unspzifsche Muskuloskeleteal backaches
 - Chr onische Neuropathiescher pain syndrome,
 - Tissues and Muskelverspannung, Biorhythmus disturbances.
 - Rule-pains, Dysmenorrhoe, among others hypertonia,
 - Because of lack at sunlight (z.b). U - boats (space travel) through too long stay in closed areas.
 - Skin-illnesses, sore-treatments.

C - contraindications:

- With pregnancy and children under 6 years should be applied the Lt hand-utensil under strict indications.
- Oversensitivity opposite-sunlight-allergy.
- Irradiation directly at the eyes and surroundings.
- with arsenic or immune Supressiva medications,
- With different paralyse-Hemiplegie in these cases passes combustion-danger.
- With acute waiter-stomach troubles, Z. B. appendicitis among others,
- With acute fever-illnesses and acute rheumatic troubles form-circle, strong headache,
- Mali-gnome, Präkanseröse, Nävus, it doesn't should

is applied.

D - piece of information **and precautions:**

- That light-thermo-therapy hand-appliance of children keeps away and unapproachable stores.
- With children under 6 years, doesn't apply.
- The application with pregnant can too premature Labor-pains is caused. During the Confinement can be applied.

E - side-effects:

- With patients with lower blood pressure, it can to Cheat attacks leads. With patients with high Blood pressure-threshold

Dear patient, dear patient!

Lead the Ltt - therapy application not through, without to have consulted previously your doctor.

You/they have the Ltt in treatment because of her/its/their troubles - light-thermo-therapy Hand-appliance gotten. Light is a most important parliamentary party of the life, already to start one,
Time without light is not conceivable the life.

The sun is source of energy the biggest highlight for us. In the winter, we get in the North-hemisphere less light, depression and pain-ailments takes through it to. Light-thermo-therapy is amicable in general very well.

Please heed the hints of the application-information and visit a doctor with prolonged troubles.

, Hypertonia, **one needs no** precautions
, **to meet. It should** with Pacemaker patients
Is **applied** carefully. With the appearance **from**
Side-effects **immediately the doctor benach right.**

F - interactions **with other** drugs:

Light-thermo-therapy hand-appliance **can with others**
Drugs is **enforced**. Her/it/them **can**
her/its/their medications **still earns**.
The **rheumatic** troubles **become fewer and**
earning of medications is **reduced**. In the
Bleed, **the fats and the cholesterol turn out right**
reduced, it hangs with the application of course
together. Light-thermo-therapy appliance **can**
sleep-promoting works.

G - type of **the application and dosage**:

- TL 12° NGES as LAFER, with ca.700 facet-mirrors, light-sources 1 cm removals per 1 cm²s 170 mWs. deliver.

Light-thermo-therapy hand-appliance **can on the day in 3**
Hours distances 1-4 **times, approximately until in each case 16 minutes**
is applied. With **very strong** pains
, Exceptions **heed** - point C, **the application is**
from light-thermo-therapy - Ltt hand-appliance **several times**
on the day, to recommend.

H-Über dosage and measures:

- **About possible combustion-dangers**
, **to avoid, asks**
we you necessarily your treating the instructions
Doctor and the instruction manual, to follow carefully.

- **With an overdose (dizziness) strong**
Headaches, black before the eyes (heavy overtiredness-appearances), stops immediately with the treatment. And your treating doctor visits.

I - storage:

- **The Ltt - hand-appliance with not-use**
exclusively in her/it original package and for children unapproachable storing.
- **The Ltt - hand-appliance before wetness and moisture protects.**
The instrument not under -10°C and over + 40°C
Area-temperature abandons. With the transportation only that original package uses.

J - manufacturers, services, offer, formula and Verpackungs - content:

One enforces manufacture and distribution through the company light-Hanseatic league in 28211 Bremen exclusively. The trade of original sharing and accessories through third-persons / companies is founded authorization of the o.a distribution on one in writing.

For repair-works, you send the hand-appliance in the ori -
ginal package with guarantee and Kassenbon at un -
sere services partners to. The sets address find you here -
to in the instruction manual on side 14.
Current offers can you in writing or telephonic
ask from our company.

With the Neukauf of the Ltt - light-thermo-therapy of hand-appliance
be in her/it original
Package:

- The Ltt - hand-appliance,
- Two pieces of energy sources,
- Small and big distance regulators per two piece,
- Service-handbook,
- Visual therapy-instruction,
- and this brochure.
- the tote bag
- Aroma - odor-stone and accessories
- Farbglässer (optional)

The Ltt - hand-appliance is under her/it authorization - and
Patent-number EP 1005 880, DP 19856002 registers.

Voucher

You have gotten a medical appliance of the company Ltd. Further Questions over the treatment like to be available we for you. With Komp - likationen you can or disturbance of the appliance itself under the towards-Ann - ten services address report.

Internet: www.lichthansa.de e-mail: [info @ lichthansa.de](mailto:info@lichthansa.de)

With this voucher, you can Patienteninfomaterial gratuitously with orders our company.

Send these at:

*Lichthansa league - med. Dr. Apul * Vionville street 7 * 28211 Bremen*

My address:

Woman / man.....Geb. Date:.....

Street:.....

PLZ / place.....

Phone, 0.....)Fax(0.....)E Mail:..... @.....